

Element HT5 Maintenance and Control Log

Month/Year:

Control lot # _____ Date opened _____
Control lot # _____ Date opened _____
Control lot # _____ Date opened _____

DAY	RUN BACKGROUND	RUN CONTROL	PROBE CLEANSER	COMMENTS	INITIALS
1	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
2	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
3	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
4	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
5	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
6	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
7	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
8	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
9	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
10	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
11	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
12	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
13	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
14	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
15	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
16	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
17	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
18	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
19	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
20	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
21	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
22	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
23	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
24	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
25	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
26	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
27	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
28	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
29	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
30	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		
31	☐ Pass ☐ Fail	☐ Pass ☐ Fail	☐		

Antech recommends performing quality control daily. Store controls at refrigerated temperature. Discard if frozen. Bring controls to room temperature and mix well prior to use. Return controls to the refrigerator immediately after use. Hematology control is stable for 14 days after opening. For further assistance, please call Technical Support Services at 800-464-3752, option 3.